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BACKGROUND

The worldwide frequency of allergies has increased significantly over the past decades. The term allergy is often used for Type I hypersensitivity reactions (immediate type reactions), whose symptoms generally occur within 30-60 minutes after contact with the allergen. The most frequent symptoms are: hay fever (rhinitis), conjunctivitis, hives (urticaria), allergic asthma and as the most dangerous manifestation anaphylaxis (the anaphylactic shock). The allergens causing Type I hypersensitivity reactions are mostly proteins derived from the natural environment e.g. plant pollen, animal hair, food, mites and insect venoms. The characteristics of Type I allergies is the involvement of allergen specific immunoglobulins (antibodies) of class E (sIgE). Hence, the detection of sIgE is an important tool of modern allergy diagnostics.

INTENDED USE

The IgE food-screening test detects specific IgE in human serum or plasma and thus adds to the diagnosis of type I allergies.

PRINCIPLE

The principle of the test is based on Enzyme Linked Immuno Sorbent Assay (ELISA) technique. Microplates were coated with food allergens of specific interest, each antigen in a separate well. Patient serum or plasma is applied to antigen coated wells. During incubation the specific IgE of the patient sample bind to the antigens (allergens). After addition of the anti-human-IgE-Enzyme-Conjugate (IgE - Horseradish - Peroxidase = IgE-HRP) a solidphase-bound antigen-IgE-Anti-IgE-HRP-complex is formed. The presence of the bound antigenantibody-complexes is detected enzymatically by application of the TMB Substrate (Tetramethylbenzidine). After stopping the enzymatic reaction with phosphoric acid (H₃PO₄) the optical density of the yellow colour, that is direct proportional to the specific IgE concentration of the patient sample, is measured spectrophotometrically at 450 nm (reference wave length 620 nm). The sIgE concentration of the patient sample is proportional to the OD. Calibrators with defined concentrations of IgE (calibrated against WHO 75/502) are assayed simultaneously with the patient samples to generate a calibration curve. Unknown IgE concentrations of the test samples are calculated from this curve.

KIT COMPONENTS

Component	10108E	10107E	10102E
Anti-IgE-Enzyme-Conjugate, ready to use	1 x 11 mL	1 x 22 mL	1 x 110 mL
Concentrated Washing Buffer	1 x 30 mL	1 x 30 mL	1 x 250 mL
TMB Substrate, ready to use	1 x 12 mL	1 x 22 mL	1 x 110 mL
Stop Solution (25% H ₃ PO ₄)	1 x 6 mL	1 x 12 mL	1 x 55 mL
Dilution Buffer	1 x 11 mL	1 x 22 mL	1 x 110 mL

MATERIAL NEEDED, BUT NOT INCLUDED IN THE KIT

1. Reference kit	12050E
Reference wells	5 x 8 wells
Calibrators	4 x 0.5 mL (17.5; 3.5; 0.7; 0.35 IU/mL)
Control high	1 x 0.5 mL
Control low	1 x 0.5 mL

2. Antigen- coated plates: REF №13035-B, №13090-A, №13074-BVR, №13075-BVR, №13076-BVR, №13085-BVR, №13040-G, №13010-N11, №13020-N12

3. LABORATORY EQUIPMENT

Pipettes: 10-100 µL, 200-1000 µL, Multipipette, pipet tips, tubes for dilution of the specimens, measuring glass, lab watch, microplate-reader, incubator, covering foil.

Optionally: microplate-washer.

SPECIMEN COLLECTION & PREPARATION

Either Serum or Plasma can be used in this test. No additives or preservatives are necessary to maintain the integrity of the specimen. Specimens should be stored at 2-8°C and assayed within 48 hours after collection. If the assay cannot be performed within 48 hours or if the specimen has to be shipped, cap the specimen and keep it frozen. Repeated freezing and thawing should be avoided. Frozen specimens should be thawed at room temperature (RT, 20-25°C) and mixed thoroughly by gentle inversion before assaying. **Samples should be tested undiluted.** The use of haemolysed or lipemic specimens is not recommended.

PREPARATION OF REAGENTS

Before starting the assay let all samples and reagents reach RT.

Dilution Buffer:	ready to use
Conjugate:	ready to use
Calibrators:	ready to use
Controls:	ready to use
TMB Substrate:	ready to use
Stop Solution:	ready to use

Concentrated Washing Buffer: Dilute concentrate Washing Buffer 1:25 in distilled water (30 mL Washing Buffer concentrate have to be diluted to a final volume of 750 mL with distilled water). The resulting Washing Buffer is stable for one week at RT.

ASSAY PROCEDURE

1. Prepare a protocol for the assay run. It is recommended to test the calibrators and controls in duplicate determination.
2. Place a reference strip and the required amount of wells for the controls in a separate frame. Reseal remaining wells with a desiccant in the aluminium bag. **Standards and controls are ready to use and have not to be diluted!** Take the appropriate number of plates for the patients from the stock.
3. Pipette 50 µL of calibrators, controls and patient samples according to the plate matrix. Cover the strips with a foil to prevent evaporation and incubate the strips for 1 h at 37 °C.
4. After the incubation time aspirate solution from the wells and wash each well 3 x with 500 µL wash solution. The washing procedure can be performed manually or by using a validated microplate-washer. Remove residual liquid by dunking the plate on a tissue.
5. Pipette 100 µL Enzym Conjugate into each well and cover the strips with a foil. Incubate 1 h at 37°C.
6. Perform the washing procedure as described in step 4.
7. Pipette 100 µL TMB Substrate into each well and cover the strips with a foil and incubate 10 min in the dark at RT.
8. Pipette 50 µL Stop Solution in the same order as the substrate to each well and mix gently. It is recommended to mix the solutions in the wells by tipping against the frame. After 5 min read the optical density using an ELISA Reader at 450 nm (Reference wavelength 620 nm).

CALCULATION OF RESULTS

It is recommended to use validated software for the calculation of the results. For manual calculation, the mean OD [D 450 nm – 620 nm] values are calculated from the calibrators. Generate a graph from the mean OD values of the four calibrators on half logarithmic paper (Abscissa: log IU IgE/mL; Ordinate: linear OD D 450 nm – 620 nm) to create a standard curve. The sIgE concentration of the patient sample is determined on the basis of this standard curve. The OD is mapped on the ordinate and the result in IU/mL can be read out on the abscissa. The standard curve and the controls should be in the acceptance range given in the Quality-Control-Certificates delivered with the kit. Otherwise, the test conditions should be verified and the test should probably be repeated.

RESULT INTERPRETATION


It is recommended that each laboratory defines and evaluates its own reference ranges.

MEASURING RANGE

Specific IgE-concentrations between 0.35 and 17.5 IU/mL can be determined.

The interpretation of the results has to be done individually for each patient and in consultation with a medical doctor or a professional nutritionist!!!

TEST SCHEME IgE

RECOMMENDATION

IgE < 0.35 IU/ml nondetectable	Foods can be eaten without problems (pay attention to the results of the IgG4-test)
0.35-0.7 IU/ml	Foods 1 to 2 times allowed per week, rotating
> 0.7 IU/ml	Avoid Foods at least 3 month (independent of IgG4-Testing)

LITERATURE

- Jarvinen KM, Chatchatee P, Bardina L, Beyer K, Sampson HA. IgE and IgG Binding Epitopes on alpha-Lactalbumin and beta-Lactoglobulin in Cow's Milk Allergy. *Int Arch Allergy Immunol* 2001; 126: 111-118
- Chatchatee P, Jarvinen KM, Bardina L, Vila L, Beyer K, Sampson HA. Identification of IgE and IgG-binding epitopes on alpha-s1-casein: Differences in patients with persistent and transient cow's milk allergy. *J Allergy Clin Immunol*. 2001; 379-383.
- Zar S, Benson MJ, Kumar D. Food-specific serum IgG4 and IgE titers to common food antigens in irritable bowel syndrome. *Am J Gastroenterol*. 2005, 100:1550-1557.

PRECAUTIONS FOR USERS

- In compliance with article 1 paragraph 2b European directive 98/79/EC the use of in-vitro diagnostic medical devices is intended to secure suitability, performance and safety of the product by the manufacturer. Therefore the test procedure, information, precautions and warnings stated in the instructions for use have to be followed strictly. The kit has only to be used as described on page 1 (intended use).
- The test must be performed according to this instruction, which contains all necessary information, precautions and warnings. The use of the test kit with analyzers and similar equipment has to be validated. Any change in design, composition of the test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes resulting in false results and other incidents. The manufacturer is not liable for any results obtained by visual analysis of patient samples.
- The kit is intended for use by trained and qualified professionals carrying out research or diagnostic activities only. Pregnant women should not perform the test.

- Laboratory equipment has to be maintained according to the manufacturer's instructions and must be tested for its correct function before use.
- For *in-vitro* diagnostic use only. Use only once. Do not use components exceeding the expiry date. Do not combine reagents of other suppliers or kit components of different lots (unless specified on page 6) with this kit.
- Do not use kit components when the package of the component is damaged. Please check all solutions prior to use for microbiological contamination. Cap vials tightly immediately after use to avoid evaporation and microbiological contamination. Do not interchange screw caps of the reagent vials.
- The kit was evaluated for use at the temperatures specified in the Testing scheme (see page 8). Higher or lower temperatures may result in values not meeting the quality control ranges.
- The washing procedure is absolutely important. Improper washing will cause erroneous results. It is recommended to use a multichannel pipette and an automated washer.
- To avoid cross-contamination and false-positive results it is recommended to perform all pipetting steps properly. Use only clean pipette tips, dispensers and lab ware.
- Test components based on human serum were tested using a CE marked method for the presence of antibodies against HIV 1 / HIV 2, Anti-HBc, and Anti-HCV as well as for hepatitis antigen HBsAg and were found to be negative. Nevertheless, material based on human serum should be handled as potentially infectious (BIOHAZARD).
- Some kit components may contain bovine serum albumin, of which according to the manufacturer no infectious potential is known. Due to the eventual occurrence of undetectable infectious agents we recommend to handle any product of animal origin as potentially infectious.
- The following safety rules should be followed with all reagents:
 - Do not get in eyes, on skin, or on clothing (P262). Do not breathe spray (P260). Pipetting should never be done by mouth, but with suitable pipetting devices.
 - IF SWALLOWED: rinse mouth. Do NOT induce vomiting (P301/330/331)
 - IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower (P303/361/353).
 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing (P303/340).
 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. (P305/351/338)
 - Don't eat, drink or smoke while performing the test. Keep away from food, feed and beverage.
 - Wear protective gloves/protective clothing/eye protection (P280). Wash hands thoroughly after handling (P264) and care for your skin.
 - Material safety data sheet is available on request.
- Stop Solution causes severe skin burns and eye damage (H314).
- TMB in high concentrations may be potentially mutagenic. Due to the low concentration of TMB in this substrate solution a mutagenic effect can be ruled out, if it is properly used.
- The preservatives (Bronidox, Thimerosal, Azid) are toxic to aquatic life, but their concentration is not hazardous to environment anymore. On disposal, flush large volumes of reagents with plenty of water. Thimerosal (WashBuf B) may cause damage to organs through prolonged or repeated exposure (H373)
- Waste containing serum must be collected in separate containers containing an appropriate disinfectant in sufficient concentration. This material has to be treated according to national biohazard and safety guidelines or regulations.
- We refer to the national regulations of medical devices regarding in-vitro diagnostic test kits.

DR. FOOKE
LABORATORIEN GMBH

Тел.: +7 (495) 799-1165
E-mail: info@fooke.ru
www.fooke.ru