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Food Intolerance Diagnostics

Food Screening IgG4	REF	10106PG	Σ	96 Determinations
Food Screening IgG4	REF	10105PG	Σ	192 Determinations
Food Screening IgG4	REF	101010PG	Σ	960 Determinations

BACKGROUND

Immunoglobulins (Ig) of the IgG4-subclass obviously play an important role in the humoral immune reaction to food antigens. Symptoms similar to those mediated by IgE have been assumed to be attributed to the presence of food specific IgG4. In so far the measurement of food specific IgG4 seems to get more importance.

INTENDED USE

The IgG4 – Food – Screening is intended to detect food antigen-specific IgG4 in human serum or plasma.

PRINCIPLE

The principle of the test is based on an Enzyme Linked Immuno Sorbent Assay (ELISA) technique. Microplates for the IgG4 Food-screening are coated with food antigens, each in a separate well. Diluted patient's serum or plasma is incubated in the antigen-coated microtiter test strips. The specific IgG4 from the patient's sample bind to the solid phase coupled antigens. After addition of anti-IgG4-Enzyme-Conjugate, a solid phase bounded antigen/IgG4/anti-IgG4-HRP-complex is formed. The amount of bound antibody corresponds to the amount of peroxidase in the respective well and is detected by incubation with the Substrate Tetramethylbenzidine (TMB) resulting in the development of a blue colour. After adding the Stop Solution the colour changes from blue to yellow. The optical density of the yellow colour is measured spectrophotometrically at 450 nm (reference wave length 620 nm). Within the measuring range the concentration of food-specific IgG4 is proportional to the intensity of the colour developed.

KIT COMPONENTS

Component	10106PG	10105PG	101010PG
Anti-IgG4-Enzyme-Conjugate, ready to use	1x 11 mL	1x 22 mL	1x 110 mL
Concentrated Washing Buffer B	1x 30 mL	1x 30 mL	1x 250 mL
TMB Substrate, ready to use	1x 12 mL	1x 22 mL	1x 110 mL
Stop Solution (25% H3PO4)	1x 6 mL	1x 12 mL	1x 55 mL
Dilution Buffer	1x 11 mL	1x 22 mL	1x 110 mL

MATERIAL NEEDED, BUT NOT INCLUDED IN THE KIT

1. Referencekit	12040PG
Reference wells	5 x 8 wells
Calibrators	7 x 0.5 mL (25, 50, 100, 250, 500, 1000, 2500 U/mL)
Control high	1x 0.5 mL
Control low	1x 0.5 mL

2. Antigen-coated plates: REF №13035-B, №13090-A, №13074-BVR, №13075-BVR, №13076-BVR, №13085-BVR, №13040-G, №13010-N11, №13020-N12

3. LABORATORY EQUIPMENT

Pipettes: 10-100 µL, 200-1000 µL, Multipipette, pipet tips, tubes for dilution of the specimens, measuring glass, lab watch, microplate-reader, incubator, covering foil.

Optionally: microplate-washer.

SPECIMEN COLLECTION & PREPARATION

Either Serum or Plasma can be used in this test. No additives or preservatives are necessary to maintain the integrity of the specimen. Specimens should be stored at 2-8°C and assayed within 48 hours after collection. If the sample has to be shipped, cap the specimen and keep it frozen. Repeated freezing and thawing should be avoided. Frozen specimens should be thawed at room temperature (RT, 20-25°C) and mixed thoroughly by gentle inversion before assaying. Samples should be diluted 1:20 (500 µL Serum + 9,5 mL Dilution Buffer). The use of haemolytic or lipemic specimens is not recommended.

PREPARATION OF REAGENTS

Before starting the assay let all samples and reagents reach RT.

Dilution Buffer:	ready to use
Conjugate:	ready to use
Calibrators:	ready to use
Controls:	to be diluted 1:20 (Dilution stable for 24 h)
TMB Substrate:	ready to use
Stop Solution	ready to use
Concentrated Washing Buffer:	Dilute 30 mL concentrated Washing Buffer with aqua bidest.). The resulting Washing Buffer is stable for one week at room temperature.

ASSAY PROCEDURE

- Follow the pipetting scheme provided with the kit.
- Dilute patient's samples and controls 1:20 with Dilution Buffer (500 µL patient's sample + 9.5 mL Dilution Buffer, 25 µL control + 475 µL Dilution Buffer). Calibrators are ready to use and have not to be diluted!
- Pipette 100 µL of the calibrators and 100 µL of the diluted controls resp. patient samples into the correspondent wells (see plate matrix).
- Cover the plate with a foil to prevent evaporation and incubate for 1 h at 37 °C.
- After the incubation time aspirate solution from all wells and wash each well 3 x with 500 µL wash solution/well. The washing procedure can be performed manually or using a validated microplate-washer. Remove residual liquid by dunking the plate on a tissue.
- Pipette 100 µL conjugate into each well and cover the plate with a foil. Incubate 1 h at 37 °C
- Perform the washing procedure as described in step 5.
- Pipette 100 µL TMB Substrate Solution into each well and cover the plate with a foil and incubate 15 min in the dark at RT.
- Pipette 50 µL Stop Solution in chronological order of the Substrat Solution into each well and mix gently. It is recommended to mix the solutions in the wells by tipping against the frame. After 5 min read the optical density using an ELISA Reader at 450 nm (Reference wavelength 620 nm) immediately.

CALCULATION OF RESULTS

It is recommended to use validated software for the calculation of the results. For manual calculation generate a graph from the OD [450 nm – 620 nm] values of the six Calibrators on half logarithmic paper (Abscissa: log U IgE/mL; Ordinate: linear OD 450 nm – 620 nm) to create a standard curve. The IgG4 concentration of the patient sample is determined on the basis of this standard curve. The OD is mapped on the Ordinate and the result in U/mL can be read out on the Abscissa. The standard curve and the controls should be in the acceptance range given in the Quality-Control-Certificates delivered with the kit. Otherwise, the test conditions should be verified and the test should probably be repeated.

Attention: Since the calibrators are not diluted it is necessary to multiply the results of the samples and the controls by 20 to obtain the specific IgG4 concentration present in the patient samples.

RESULT INTERPRETATION

It is recommended that each laboratory defines and evaluates its own reference ranges.


MEASURING RANGE

This measuring range detects IgG4 concentrations between 500 U/mL and 50.000 U/mL.

RECOMMENDATION

< 2000 U/mL IgG4	Foods can be eaten without limitation Pay attention to possible IgE results
> 2000 < 5000 U/mL IgG4	Foods allowed 1 to 2 times per week, rotating Pay attention to possible IgE results
> 5000 U/mL IgG4	Avoid Foods at least 3 month Pay attention to possible IgE results

1 U/mL ≈ 1 ng/mL

TEST SCHEME IgG4

The interpretation of the results has to be done individually for each patient and in consultation with a medical doctor or a professional nutritionist!!!

PRECISION

Variability and Reproducibility
Intra-Assay-Variability

Specimen	Mean [IU/mL]	CV (%)
1 (n=3)	>3000	<11
2 (n=3)	>3000	<21

LINEARITY

Four randomly selected sera show a linear behaviour (< ±20%) in three consecutive dilution steps. Based on the heterogeneity of human serum or plasma samples varying results can not be excluded.

LITERATURE

- Jarvinen KM, Chatchatee P, Bardina L, Beyer K, Sampson HA: IgE and IgG Binding Epitopes on alpha-Lactalbumin and beta-Lactoglobulin in Cow s Milk Allergy. *Int Arch Allergy Immunol* 2001, 126:111-118.
- Chatchatee P, Jarvinen KM, Bardina L, Vila L, Beyer K, Sampson HA: Identification of IgE- and IgG-binding epitopes on alpha(s1)-casein: Differences in patients with persistent and transient cow s milk allergy. *J Allergy Clin Immunol* 2001, 107:379-383.
- Kruszewski J, Raczka A, Klos M, Wiktor-Jedrzejczak W: High Serum Levels of Allergen Specific IgG-4 (asIgG-4) for Common Food Allergens in Healthy Blood Donors. *Archivum Immunologiae et Therapiae Experimentalis* 1994, 42:259-261.
- Zar S, Mincher L, Benson MJ, Kumar D: Food-specific IgG4 antibody-guided exclusion diet improves symptoms and rectal compliance in irritable bowel syndrome. *Scand J Gastroenterol* 2005, 40:800-807.
- Zar S, Benson MJ, Kumar D: Food-specific serum IgG4 and IgE titers to common food antigens in irritable bowel syndrome. *Am J Gastroenterol* 2005, 100:1550-1557.
- Noh G, Ahn HS, Cho NY, Lee S, Oh JW: The clinical significance of food specific IgE/IgG4 in food specific atopic dermatitis. *Pediatr Allergy Immunol* 2007, 18: 63-70

PRECAUTIONS FOR USERS

- In compliance with article 1 paragraph 2b European directive 98/79/EC the use of in-vitro diagnostic medical devices is intended to secure suitability, performance and safety of the product by the manufacturer. Therefore the test procedure, information, precautions and warnings stated in the instructions for use have to be followed strictly. The kit has only to be used as described on page 1 (intended use).
- The test must be performed according to this instruction, which contains all necessary information, precautions and warnings. The use of the test kit with analyzers and similar equipment has to be validated. Any change in design, composition of the test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes resulting in false results and other incidents. The manufacturer is not liable for any results obtained by visual analysis of patient samples.
- The kit is intended for use by trained and qualified professionals carrying out research or diagnostic activities only. Pregnant women should not perform the test.
- Laboratory equipment has to be maintained according to the manufacturer's instructions and must be tested for its correct function before use.
- For in-vitro diagnostic use only. Use only once. Do not use components exceeding the expiry date. Do not combine reagents of other suppliers or kit components of different lots (unless specified on page 6) with this kit.
- Do not use kit components when the package of the component is damaged. Please check all solutions prior to use for microbiological contamination. Cap vials tightly immediately after use to avoid evaporation and microbiological contamination. Do not interchange screw caps of the reagent vials.
- The kit was evaluated for use at the temperatures specified in the Testing scheme (see page 8). Higher or lower temperatures may result in values not meeting the quality control ranges.
- The washing procedure is absolutely important. Improper washing will cause erroneous results. It is recommended to use a multichannel pipette and an automated washer.
- To avoid cross-contamination and false-positive results it is recommended to perform all pipetting steps properly. Use only clean pipette tips, dispensers and lab ware.
- Test components based on human serum were tested using a CE marked method for the presence of antibodies against HIV 1 / HIV 2, Anti-HBc, and Anti-HCV as well as for hepatitis antigen HBsAg and were found to be negative. Nevertheless, material based on human serum should be handled as potentially infectious (BIOHAZARD).
- Some kit components may contain bovine serum albumin, of which according to the manufacturer no infectious potential is known. Due to the eventual occurrence of undetectable infectious agents we recommend to handle any product of animal origin as potentially infectious.
- The following safety rules should be followed with all reagents:
 - Do not get in eyes, on skin, or on clothing (P262). Do not breathe spray (P260). Pipetting should never be done by mouth, but with suitable pipetting devices.
 - IF SWALLOWED: rinse mouth. Do NOT induce vomiting (P301/330/331)
 - IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower (P303/361/353).
 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing (P303/340).
 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. (P305/351/338)
 - Don't eat, drink or smoke while performing the test. Keep away from food, feed and beverage.
 - Wear protective gloves/protective clothing/eye protection (P280). Wash hands thoroughly after handling (P264) and care for your skin.
 - Material safety data sheet is available on request.
- Stop Solution causes severe skin burns and eye damage (H314).
- TMB in high concentrations may be potentially mutagenic. Due to the low concentration of TMB in this substrate solution a mutagenic effect can be ruled out, if it is properly used.
- The preservatives (Bronidox, Thimerosal, Azid) are toxic to aquatic life, but their concentration is not hazardous to environment anymore. On disposal, flush large volumes of reagents with plenty of water. Thimerosal (WashBuf B) may cause damage to organs through prolonged or repeated exposure(H373)
- Waste containing serum must be collected in separate containers containing an appropriate disinfectant in sufficient concentration. This material has to be treated according to national biohazard and safety guidelines or regulations.
- We refer to the national regulations of medical devices regarding in-vitro diagnostic test kits.